Zählwaage / Counting scale / Balance compteuse

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Supplemental directives

Read this document carefully and completely before using the device and observe the safety precautions contained therein. The manufacturer is not liable for damage resulting from disregarding the instructions or improper use.

Table of contents

1	Introduction.....	113	5.2.3	Removing the transport safety lock....	132
1.1	Validity.....	113	5.2.4	Levelling the device.....	133
1.2	Document purpose.....	113	5.2.5	Connecting wired peripheral devices..	134
1.3	Other versions of the document.....	114	5.2.6	Mains operation.....	134
1.4	Presentation conventions.....	114	5.2.7	Operation with rechargeable battery pack.....	136
2	General safety instructions.....	116	5.3	Commissioning the device.....	138
2.1	Meaning of the warnings.....	116	6	Settings.....	139
2.2	Safe handling of packaging material...	117	6.1	Opening and using the setup menu....	139
2.3	Safe handling of electrical devices.....	117	6.2	Setup menu overview.....	140
2.4	Notes on IT security.....	120	6.3	Adjustment settings.....	141
2.5	Intended use.....	120	6.3.1	Adjusting procedure with an external adjustment weight.....	142
2.6	Improper use.....	120	6.3.2	Adjusting procedure with a user- defined adjustment weight.....	143
2.7	Obligations of the operating company	120	6.4	Setting the acoustic signal.....	144
2.8	Required user qualifications.....	121	6.5	Setting the automatic switch-off func- tion.....	144
2.9	Requirements for the environment of use.....	121	6.6	Setting the automatic switch-on func- tion.....	145
3	Description of the device.....	123	6.7	Assigning of key functions.....	145
3.1	General product description.....	123	6.8	Adjusting the backlight.....	146
3.2	Overview of the device.....	123	6.9	Setting the taring range.....	146
3.3	Overview of the keys.....	124	6.10	Setting the zero range (Zero Limit)....	146
3.4	Overview of the display.....	126	6.11	Setting automatic zero tracking.....	147
4	Unpacking, checking and trans- porting the device.....	128	6.12	Setting the display speed / filter.....	147
4.1	Unpacking and checking the device...	128	6.13	Setting the availability of display units	147
4.2	Scope of delivery.....	129	6.14	Setting the availability of applications	147
4.3	Transporting the device.....	129	6.15	Restoring factory settings.....	148
5	Assembly, installation and commis- sioning.....	131	7	Interfaces and data output.....	149
5.1	Selecting the place of use.....	131	7.1	Optional software.....	149
5.2	Installing the device.....	131	7.2	Description of the interfaces.....	149
5.2.1	Safety during installation.....	131			
5.2.2	Preparing for installation.....	131			

Table of contents

7.3	Setting up interfaces.....	149	11 Tolerance weighing application.....	181
7.3.1	Connecting wired peripheral devices..	149	11.1	Overview of the "Tolerance weighing" application menu.....
7.3.2	Setting interface parameters.....	150	11.2	Weighing the target weight.....
7.4	Configuring data output.....	151	11.3	Weight checking (limit values).....
7.4.1	Overview of output settings.....	151	12 Other features.....	186
7.4.2	Setting output parameters.....	154	12.1	Applying the totalizing function.....
7.4.3	Print settings.....	157	13 Care, maintenance and servicing.....	187
7.5	KCP - KERN Communications Protocol (KERN interface protocol).....	160	13.1	Cleaning.....
8 Basic operation.....	161		13.1.1	Safety during cleaning.....
8.1	Switching on the device.....	161	13.1.2	Recommended cleaning.....
8.2	Switching off the device.....	162	13.2	Servicing.....
8.3	Zero setting.....	162	13.2.1	Safety during servicing.....
8.4	Taring.....	163	13.2.2	Maintenance schedule.....
8.4.1	Tare.....	163	13.3	Control of inspection, measuring and test equipment.....
8.4.2	Saving a known tare weight (pre-tare)	164	14 Decommissioning and storage.....	190
8.5	Selecting an application.....	165	14.1	Taking the device out of operation.....
9 Weighing application.....	167		14.2	Storing the device.....
9.1	Overview of the "Weighing" application menu.....	167	15 Warranty and goods return.....	191
9.2	Simple weighing.....	167	15.1	Warranty.....
9.3	Underfloor weighing.....	167	15.2	Goods return.....
9.4	Hold function.....	168	16 Disposal.....	192
9.5	Percentage weighing.....	169	17 Troubleshooting.....	194
9.6	Weighing with a multiplication factor..	170	17.1	Error messages.....
10 Counting application.....	172		17.2	Troubleshooting table.....
10.1	Overview of the "Counting" application menu.....	172	18 Product specification.....	197
10.2	Setting the reference quantity.....	174	18.1	Technical data.....
10.3	Simple piece counting.....	176	18.2	Declaration of Conformity.....
10.4	Weighing the target number of items..	177	19 Contact details.....	202
10.5	Quantity checking (limit values).....	179	20 Index.....	203

1 Introduction

1.1 Validity

This document is an integral part of the product. It is valid for the product versions listed below.

You can recognise which product version your product is on the type plate. The following illustration is an example. The real type plate may differ from the illustration. The model number is shown in <Model> and the type item no. is shown in <Type>.


Fig. 1: Example of the type plate (real type plate may be different)

Model	Item no. / Type
CDS 4K0.02	TCDS 4K0.02-A
CDS 15K0.05	TCDS 5K0.05-A
CDS 16K0.1	TCDS 16K0.1-A
CDS 30K0.1	TCDS 30K0.1-A
CDS 30K0.1L	TCDS 30K0.1L-B
CDS 36K0.2L	TCDS 36K0.2L-B
CDS 60K0.2	TCDS 60K0.2-B

1.2 Document purpose

This document provides the user with descriptions for the safe and efficient use of the product. The document must be carefully read and understood by all users before the product is used.

1.3 Other versions of the document

The original version of this document was created in German. We reserve the right to make subsequent changes to the product documentation. The current version and other language versions of the document can be found on our website www.kern-sohn.com.

You can find more interesting information about our products in the download area:

- Technical data sheets
- Product brochures and catalogues
- Information on the KCP protocol and software

1.4 Presentation conventions

References in the text


Element	Meaning
Example: [TARE]	Labelling of keys
Example: <5E6.0P>	Text on the scale display
*	Asterisk * indicates default settings

Instructions

Instructions are presented as step-by-step guidance.

1.  This is step 1.
2.  This is step 2.
 - ➡ This is the result of step 2.

Scale operation

Element	Meaning
	left = action, right = result (the scales and displays shown may differ from this example)
	Navigation keys for input on the device (the keys shown may differ from those in this example)

Other representations

Here you will find additional information or tips.

2 General safety instructions

Read this document thoroughly and carefully before using the device and observe all safety precautions contained herein to avoid bodily injury and damage to property. Keep this document in a place that is always easily accessible so that it is quickly available when needed.

2.1 Meaning of the warnings

Warnings warn you of possible safety risks during certain activities. Read the warnings in full and follow the instructions. Failure to do so may result in injury, damage to property, environmental damage, malfunctions or incorrect results.

Meaning of the signal words

Signal word	Meaning
DANGER	This combination of symbol and signal word indicates an immediate dangerous situation that will result in death or serious injury if it is not avoided.
WARNING	This combination of symbol and signal word indicates a possible dangerous situation that can result in death or serious injury if it is not avoided.
CAUTION	This combination of symbol and signal word indicates a possible dangerous situation that can result in minor injury if it is not avoided.
NOTICE	This combination of symbol and signal word indicates a possible dangerous situation that can result in material and environmental damage if it is not avoided.
ENVIRONMENT	This combination of symbol and signal word indicates a possible dangerous situation that can result in material and environmental damage if it is not avoided.

Meaning of the warning signs

Warning signs	Type of danger
	Warning – danger zone.
	Warning – dangers due to batteries.
	Warning – high-voltage.

2.2 Safe handling of packaging material

Always keep packaging material away from children and animals to avoid the danger of suffocation.

2.3 Safe handling of electrical devices

Safe handling of power supply units and power supplies

Improper handling of mains-powered devices or the use of damaged or unsuitable power supply units or mains cables can lead to life-threatening injuries. Therefore, please observe the following instructions:

- If the electrical equipment is defective (e.g. damage to the power supply unit or mains cable, deformation, smoke, etc.), the device must not be used. If you notice any damage, disconnect the device from the power supply immediately and contact your dealer or the technical service of KERN & SOHN.
- Only connect the device to a properly installed and undamaged power supply. The information on the device (type plate) and the mains voltage at the place of use must be identical.
- Only use original power supply units and connecting cables. If you need a replacement, please contact your dealer or the technical service of KERN & SOHN.

- Any work on the electrical equipment of the device may only be carried out by the technical service of KERN & SOHN.
- The mains cable must not be pinched or kinked and must be laid in such a way that it does not pose a tripping hazard.
- In case of an emergency, it must be possible to disconnect the device from the power supply without any problems.

Safe handling of batteries and rechargeable batteries

Batteries and rechargeable batteries may catch fire, emit toxic fumes, release corrosive liquids or explode if not used correctly. Therefore, observe the instructions in the following section.

Damage

- Do not use damaged or deformed batteries and rechargeable batteries.
- Protect batteries and rechargeable batteries from mechanical damage, heat, cold, water and chemicals.
- Protect batteries and rechargeable batteries from short circuits. Never touch the terminals with metal objects.
- When inserting batteries and rechargeable batteries, ensure that the polarity (positive and negative terminals) is correct and that you are using suitable battery types. Only use types recommended by the manufacturer.
- Batteries and rechargeable batteries must not be modified.

Charging batteries

- Chargers must not be modified.
- Only use the original charger.
- Observe the charging time specified in the technical data and avoid overcharging.
- Do not charge batteries that are deformed or show other signs of damage.
- In the event of heat build-up, deformation, smoke or burnt smell, disconnect the charger immediately and remove it from service.

Transport

- When transporting batteries or rechargeable batteries, the applicable national and international regulations must be observed. Ensure that batteries or rechargeable batteries are protected against short circuits and damage. Do not ship defective batteries or devices containing defective batteries.

Storage

- Charge the battery if you will not be using it for a long time.
- Remove the batteries or rechargeable batteries if the device is not to be used for a longer period of time.
- Protect batteries and rechargeable batteries from short circuits (e.g. by covering the terminals) and store them separately by type. Preferably use a fireproof container for storage.
- Store batteries and rechargeable batteries at room temperature in a dry place. Avoid extreme cold or extreme heat in the storage environment. Also observe the recommended ambient temperature stated in the technical data.

Battery fluid

- Avoid contact between leaking fluid and skin, eyes or clothing.
- Wear protective clothing/equipment when removing a defective battery and ensure the room is well ventilated.
- Clean skin or clothing that has come into contact with battery fluid thoroughly with soap suds and then rinse with clean water.
- In case of contact with eyes, rinse immediately with clean water and then get medical help.

Disposal

- Do not dispose of batteries and rechargeable batteries in household waste, but at collection points or recycling centres.
- Ensure that the poles are not short-circuited during disposal.
- Further information on disposal: ➔ *'Disposal of batteries and rechargeable batteries' on page 192*

2.4 Notes on IT security

This device is equipped with an interface for data transmission. Unauthorised persons may be able to gain access to the data being exchanged. It is the operator's responsibility to identify the potential IT risks at the place of use and to take appropriate countermeasures where necessary.

2.5 Intended use

Intended use of the scale

This scale is used to determine the weight of items to be weighed. This scale is a "non-automatic weighing instrument". This means that the material to be weighed must be placed manually, with care and as centred as possible on the weighing plate. Once a stable weight value is displayed, the result can be read.

2.6 Improper use

Improper use of the scale

- Our scales with strain gauges are not intended for use in dynamic weighing processes.
However, the scales can also be used for dosing applications after checking the individual area of application and the accuracy requirements.
You can find out whether your scale is equipped with strain gauge technology by checking the technical data.
➔ *Chapter 18.1 'Technical data' on page 197*
- To prevent damage, avoid impact loadings and do not exceed the maximum load limit (Max.) of the scale.
- Any modification to the design of the scale will invalidate the warranty and may lead to incorrect weighing results and safety-related defects. Such modifications are therefore prohibited.
- The scale may only be used in accordance with the specifications described in this document.
- The scale is not designed for weighing people. It is not suitable for medical purposes.

2.7 Obligations of the operating company

It is the responsibility of the operating company to comply with the specifications described in this document as well as with legal and official regulations at the place of use. Any failure to comply with the safety precautions and instruc-

tions described in this document is considered improper use. KERN & SOHN is not liable for any damage resulting from improper use.

2.8 Required user qualifications

- Assembly, commissioning, servicing and maintenance of the device may only be carried out by qualified personnel. A qualified person is a person with appropriate professional training or qualifications for the activities described. They have in-depth knowledge of the tasks and the relevant safety standards.
- Only trained users may operate, handle and maintain the device. A trained user is a person who has been trained in the use of the device and in the corresponding safety instructions.

2.9 Requirements for the environment of use

In addition to the environmental information in the technical data, observe the following requirements for the environment of use:

- The device must not be operated in environments where there is a risk of corrosion.
- The device must be protected from excessive humidity, steam, spray and dust.
- The device must be protected from excessive moisture caused by condensation. Condensation can occur if the device is moved from a cold environment to a much warmer environment. If such a change of environment cannot be avoided, allow the device to acclimatise for approx. 2 hours before starting it up.
- The device must not be used in environments with aggressive chemicals, gases or liquids.
- The device must not be used in potentially explosive atmospheres or environments where explosives are present.
- The device must not be used in environments subject to severe impacts or vibrations.
- The device must not be exposed to extreme heat or strong temperature fluctuations (e.g. by placing the device next to a heater or in direct sunlight).
- Protect the device from direct draughts.

- Select a location with a stable power supply.
- When choosing the place of use, consider the surrounding environment. Electromagnetic fields and electrostatic charges can affect the weighing result or cause damage.

3 Description of the device

3.1 General product description

This scale is a counting scale. Counting scales are particularly suitable for use in production facilities to calculate the number of units produced based on their weight.

3.2 Overview of the device


- 1 Dimensions of the weighing plate
- 2 Mains adapter connection
- 3 KUP connection (KERN Universal Port)
- 4 Display

- 5 Keys
- 6 Battery compartment
- 7 Level indicator

Description of the device

		
CDS 4K0.02	CDS 15K0.05 CDS 16K0.1 CDS 30K0.1	CDS 30K0.1L CDS 36K0.2L CDS 60K0.2

3.3 Overview of the keys

Key	Designation in the text	Short keystroke	Long keystroke
1	[ON/OFF] or [←]	<ul style="list-style-type: none"> ■ Switch on the scale ■ Switch backlight on/off (only on devices with backlighting) ■ Move up one menu level / Exit menu ■ Previous digit position 	<ul style="list-style-type: none"> ■ Switch off the scale

Key	Designation in the text	Short keystroke	Long keystroke
2	[TARE] or [→]	<ul style="list-style-type: none"> ■ Tare ■ Select menu item and move to the next menu level / Confirm selection ■ Confirm digit and move to the next digit position 	<ul style="list-style-type: none"> ■ Open the application menu
3	[5x]	Reference quantity "5"	
4	[10x/Ref n]	Reference quantity "10" (short keystroke) or Select a free reference quantity (long keystroke)	
5	[20x]	Reference quantity "20"	
6	[CHANGE] or [↓]	<ul style="list-style-type: none"> ■ Next menu item ■ Decrease digit <p>Application: Weighing</p> <ul style="list-style-type: none"> ■ Switch between weighing units (when pressed for the first time: set the secondary weighing unit to be displayed) <p>Application: Counting</p> <ul style="list-style-type: none"> ■ Switch between number of pieces and weight (when pressed for the first time: set the secondary weighing unit to be displayed. Note: If pcs is selected, a reference must then be established.) 	<p>Application: Weighing</p> <ul style="list-style-type: none"> ■ Display gross weight ■ Display tare weight <p>Application: Counting</p> <ul style="list-style-type: none"> ■ Display gross weight or gross quantity ■ Display tare weight or tare quantity <p>Application: Tolerance weighing</p> <ul style="list-style-type: none"> ■ Switch the display between gross, net and tare weight

Key	Designation in the text	Short keystroke	Long keystroke
		Application: Tolerance weighing <ul style="list-style-type: none"> ■ Switch between weighing units (when pressed for the first time: set the secondary weighing unit to be displayed) 	
7	[PRINT] or [↑]	<ul style="list-style-type: none"> ■ Send data to a connected peripheral device ■ Previous menu item ■ Increase digit 	

3.4 Overview of the display

Display / Icon	Designation	Description
	Indicator: Stability	Displayed when the weight value is stable
>0<	Indicator: Zero	Displayed while the scale is at zero
	Indicator: Negative	Displayed when the weight reading is negative
	"HI" tolerance mark	Upper tolerance limit exceeded
	"OK" tolerance mark	Weight reading within the tolerance limits
	"LO" tolerance mark	Value below the lower tolerance limit
PTARE	Indicator: Pre-tare	Displayed when pre-tare is active
HOLD	Indicator: Hold	Displayed when weighing values are held
NET	Indicator: Net	The displayed value is the net weight (shown after taring)

Display / Icon	Designation	Description
G	Indicator: Gross weight	The displayed value is the gross weight (shown after zero setting)
AP	Indicator: Auto-print	Displayed when the auto-print function is active
Σ	Indicator: Total	Displayed when the totalizing function is used
	Indicator: Data transmission	Displayed during data transfer
	Indicator: Wi-Fi	Displayed when connected to Wi-Fi
	Character on the main display	The main display shows the results or messages generated by the scale
<p>pcslbFFA %ctg/mol kgmgn dwtoztli_h_s</p>	Units	Active unit displayed on the right-hand side of the display
	Battery charging indicator	Displayed when running on battery power

4 Unpacking, checking and transporting the device

4.1 Unpacking and checking the device

WARNING

Injuries due to heavy weight of the device

Lifting heavy loads can lead to injuries and cause permanent health damage (e.g. damage to the spine)

- Only remove the device from the packaging with the help of another person.
- Use a suitable lifting device to remove or transport the device. Only use lifting equipment if you are appropriately qualified and trained. If you do not have the necessary qualifications, ask trained persons to help you. Secure the device properly.

WARNING

Danger of suffocation due to packaging material

Unattended packaging material can become a suffocation hazard for children or animals.

- Do not leave packaging material unattended, especially if children or animals are around.

- *Save the original packaging and its components for transport, storage or a possible return.*
- *Contact your dealer if you notice any missing parts or damage. Information on warranty and goods return can be found in ➔ Chapter 15.1 'Warranty' on page 191.*

1. ➔ Check the packaging for visible external damage that may have occurred during shipment.

2. ➤ Unpack all parts of the delivery (➔ *Chapter 4.2 'Scope of delivery' on page 129*) and remove the packaging materials.
3. ➤ Check that all parts of the delivery are complete and undamaged.
 - ➔ You can now transport the device to the place of use.
➔ *Chapter 4.3 'Transporting the device' on page 129*

4.2 Scope of delivery

- Scale
- Power supply unit
- Power plug adapter set
- Protective working cover
- Operating instructions

4.3 Transporting the device

WARNING

Injuries due to heavy weight of the device

Lifting heavy loads can lead to injuries and cause permanent health damage (e.g. damage to the spine)

- Only remove the device from the packaging with the help of another person.
- Use a suitable lifting device to remove or transport the device. Only use lifting equipment if you are appropriately qualified and trained. If you do not have the necessary qualifications, ask trained persons to help you. Secure the device properly.

NOTICE

Damage due to improper transport

Incorrect transport can damage the device.

General information on transporting scales

- Remove all loads from the weighing plate before transporting the scale.
- Remove all loose components from the device before transporting it.
- Re-attach the transport safety lock before transporting the device. This applies in particular to transport in vehicles.
- Use the original packaging to transport the device over long distances.
- When transporting batteries or rechargeable batteries, the applicable national and international regulations must be observed. Ensure that batteries or rechargeable batteries are protected against short circuits and damage. Do not ship defective batteries or devices containing defective batteries.
- When transporting the scale, hold it by the housing and not by the weighing plate or the weighing plate support.
- Use suitable lifting and transport equipment to move the scale. Secure the scale before transport to prevent it from falling.

5 Assembly, installation and commissioning

5.1 Selecting the place of use

The scales are designed to achieve reliable results under normal operating conditions. You can work accurately and quickly if you choose the right location for using your scale.

Please note when choosing the location:

- ➔ *Chapter 18.1 'Technical data' on page 197*
- ➔ *Chapter 2.9 'Requirements for the environment of use' on page 121*

5.2 Installing the device

5.2.1 Safety during installation

DANGER

Electric shock due to live parts

Contact with live parts can lead to electric shock and burns.

- Installation, cleaning, maintenance and servicing work may only be carried out by qualified personnel.
- Disconnect live parts from the power supply.
- Only carry out installation, cleaning, servicing and maintenance work on devices that are switched off and disconnected from the power supply, unless otherwise specified in the relevant instructions.
- Follow the instructions in this document.

5.2.2 Preparing for installation

1. ➔ Take the device to the place of use (➔ *Chapter 4.3 'Transporting the device' on page 129*, ➔ *Chapter 2.9 'Requirements for the environment of use' on page 121*)
2. ➔ Ensure that all parts included in the scope of delivery, any tools that may be required, and the relevant documentation are at hand.



5.2.3 Removing the transport safety lock

NOTICE

Do not remove any screws other than those of the transport safety lock, as this may damage the scale.

- Remove the transport safety lock as shown in the table below.

<p>Models</p> <ul style="list-style-type: none"> ■ CDS 4K0.02 	<p>Models</p> <ul style="list-style-type: none"> ■ CDS 30K0.1L ■ CDS 36K0.2L ■ CDS 60K0.2
<p>Transport lock with one screw: The transport safety lock is located on the bottom side of the platform and is marked with an arrow.</p> 	<p>Transport lock with four screws: The four screws 1 are located underneath the weighing plate.</p>  <p>The corner load stops 2 must not be removed.</p>


5.2.4 Levelling the device

- Select a suitable place of use for the scale.
➔ Chapter 2.9 'Requirements for the environment of use' on page 121
- Check the levelling regularly.
- Level the scale after moving it to a new location.

For accurate measurements, it is necessary for the scale to be levelled at its place of use.

- 1.** ➔ Position the scale at the place of use.

- 2.** ➔ Adjust the levelling feet **1** until the air bubble of the level indicator **2** is in the prescribed circle.

➔ The scale is now ready for use.

5.2.5 Connecting wired peripheral devices

*Please use only KERN accessories and printers.
These are perfectly suited to your scale.*


Unplug the scale from the mains before connecting any peripheral devices (e.g. a printer or a PC).

For information on configuring the interfaces, see ➔ [Chapter 7 'Interfaces and data output'](#) on page 149.

5.2.6 Mains operation

Fitting the power plug adapter

A power plug adapter set is included with the device so that it can be used flexibly in different locations with different power plug standards. The installation of the mains plug adapter is described below.

1. ➔ Press the release tab 1.
2. ➔ Slide the power plug adapter 2 into the power supply unit using the guide until it clicks into place.
 - ➔ The power supply unit is now ready for use.
 - ➔ To remove the power plug adapter, press the release tab and proceed in reverse order.

Connecting the power supply to the mains

WARNING

Serious injuries due to damaged power supply units

Damaged power supply units can lead to short circuits, fires and serious injuries.

- Before each use, check the power supply unit and mains cable for visible damage. Do not use the power supply unit if it is damaged.

WARNING

Injuries due to incorrect handling of power supply units

Incorrect handling of power supply units can lead to short circuits, fires and serious injuries.

- Observe the safety precautions for handling power supply units. ➔ *'Safe handling of power supply units and power supplies' on page 117*
- Avoid moisture at the mains connection, power supply unit and socket.
- Only connect the device to a properly installed socket.
- Make sure that the power supply unit cable is never kinked or pinched.
- Ensure that the power supply unit is accessible at all times.
- Ensure that the mains cable does not pose a tripping hazard.
- Only connect the power supply unit to the mains if the information on the device (sticker) and the mains voltage at the place of use are identical.

1. ➔ Plug the mains cable connection into the mains connection of the device (for the mains connection, see ➔ *Chapter 3.2 'Overview of the device' on page 123*).
2. ➔ Insert the power plug into the socket.
 - ➔ The device can now be put into operation. ➔ *Chapter 5.3 'Commissioning the device' on page 138*

5.2.7 Operation with rechargeable battery pack

WARNING

Injuries due to incorrect handling of batteries or rechargeable batteries

Incorrect handling of batteries can lead to short circuits, fires, leaks and injuries.

- Never use batteries or rechargeable batteries that are damaged or deformed.
- Never touch the terminals of batteries or rechargeable batteries with metal objects, as this can lead to short circuits.
- When inserting batteries or rechargeable batteries, ensure that the polarity is correct (positive and negative polarity).
- Do not use different types of batteries or old and new batteries together if several batteries are required for operation. This can lead to damage.
- Only replace batteries or rechargeable batteries with types recommended by the manufacturer.
- Remove the batteries from the device if you wish to store them for a longer period of time. This protects the device from damage caused by any leakage of battery fluid.

The rechargeable battery is available as an optional extra. For further information, please visit our website at ➔ www.kern-sohn.com.

Insert the rechargeable battery

1. ➤ Switch off the device and disconnect it from the mains.
2. ➤ Open the battery compartment cover on the back side of the display device.
3. ➤ Connect the rechargeable battery to the pin connector and place it in the battery compartment.
4. ➤ Close the battery compartment.

Recharge the battery

WARNING

Serious injuries due to damaged power supply units/chargers

Damaged power supply units/chargers can lead to short circuits, fires and serious injuries.

- Before each use, check the power supply unit/charger and mains cable for any visible damage. Do not use the power supply unit/charger if it is damaged.

WARNING

Injuries due to incorrect handling of power supply units/chargers

Incorrect handling of power supply units/chargers can lead to short circuits, fires and serious injuries.

- Observe the safety precautions for handling power supply units/chargers. ➔ *'Safe handling of power supply units and power supplies' on page 117*
- Avoid moisture at the mains connection, power supply unit/charger and at the socket.
- Only connect the device to a properly installed socket.
- Make sure that the power supply unit/charger cable is never kinked or pinched.
- Make sure that the power supply unit/charger is accessible at all times while the battery is charging.
- Ensure that the mains cable does not pose a tripping hazard.
- Only connect the power supply unit/charger to the mains if the specifications on the device (sticker) match the mains voltage at the place of use.

WARNING

Injuries caused by improper charging of rechargeable batteries

If batteries are not charged correctly, they may catch fire or explode.

- Chargers must not be modified.
- Only use the original charger.
- Observe the charging time specified in the technical data and avoid overcharging.
- Do not charge batteries that are deformed or show other signs of damage.
- In the event of heat build-up, deformation, smoke or burnt smell, disconnect the charger immediately and remove it from service.
- Please observe the general safety instructions regarding batteries and rechargeable batteries. ➔ *'Safe handling of batteries and rechargeable batteries' on page 118*

The battery should be fully charged before first use. You can find the charging time in the technical data ➔ *Chapter 18.1 'Technical data' on page 197.*

- To charge the battery, connect the power supply or charger to the scale.

5.3 Commissioning the device

Measures before commissioning

In order to obtain accurate weighing results with electronic scales, the scale must have reached its operating temperature. During this acclimatisation period (➔ *Chapter 18.1 'Technical data' on page 197*), the scale must be connected to the power supply.

Commissioning procedure


Complete all necessary installation work and preparations before commissioning.

1. ► Switch on the device. ➔ *Further information on page 162*
2. ► An adjustment may be required at the installation location at the beginning. ➔ *Chapter 6.3 'Adjustment settings' on page 141*

6 Settings

6.1 Opening and using the setup menu

Open the setup menu

- ▶ In operating mode, press and hold **[TARE]** and **[ON/OFF]** simultaneously.

➔ <SETUP> → <CAL> is displayed.

➔ The scale is in setup menu mode.

Menu navigation

Key	Function
[↑]	Previous menu item
[↓]	Next menu item
[→]	Select menu item and move to the next menu level / Confirm selection
[←]	Move up one menu level / Exit menu

Numerical input

Key	Function
[↑]	Increase digit
[↓]	Decrease digit
[→]	Confirm digit and move to the next digit position
[←]	Previous digit position

6.2 Setup menu overview

Overview of the main menu

The setup menu consists of several levels. The settings apply globally to all applications.

Menu level	Description	Further information
<CAL>	Adjusting	➔ Chapter 6.3 'Adjustment settings' on page 141
<CONF>	Interface parameters	➔ Chapter 7.3.2 'Setting interface parameters' on page 150
<Print>	Data output and printing	➔ Chapter 7.4.1 'Overview of output settings' on page 151
<KEYTON>	Keypad sound	➔ Chapter 6.4 'Setting the acoustic signal' on page 144
<AutoFF>	Automatic shut-off function	➔ Chapter 6.5 'Setting the automatic switch-off function' on page 144
<CAr>	Automatic switch-on function	➔ Chapter 6.6 'Setting the automatic switch-on function' on page 145
<button>	Assignment of keys	➔ Chapter 6.7 'Assigning of key functions' on page 145
<BLIGHT>	Backlighting	➔ Chapter 6.8 'Adjusting the backlight' on page 146
<TAREr>	Taring range	➔ Chapter 6.9 'Setting the taring range' on page 146
<ZLIT>	Zero-setting range	➔ Chapter 6.10 'Setting the zero range (Zero Limit)' on page 146
<ZTrAcF>	Automatic zero tracking	➔ Chapter 6.11 'Setting automatic zero tracking' on page 147
<FILTEr>	Speed of display	➔ Chapter 6.12 'Setting the display speed / filter' on page 147
<unIt>	Available units	➔ Chapter 6.13 'Setting the availability of display units' on page 147

Menu level	Description	Further information
<MODE>	Available applications	➔ Chapter 6.14 'Setting the availability of applications' on page 147
<RESET>	Restore factory settings (except date and time)	➔ Chapter 6.15 'Restoring factory settings' on page 148

6.3 Adjustment settings

The acceleration due to gravity varies at different locations on earth. For exact weighing results, the scale must be adjusted to the prevailing acceleration due to gravity at the location where it is installed. The scale must have reached its operating temperature before it can be adjusted.

Adjustment of the scale is required in the following cases:

- Before initial commissioning, if the scale has not already been adjusted to the installation location at the factory.
- After every change of location.
- After a significant change in environmental conditions (e.g. temperature or humidity).
- If the scale has been disconnected from the power supply.


To ensure accurate measurements, it is also recommended that the scale is adjusted regularly.

**NOTICE****Incorrect weighing results**

- Avoid unstable environmental conditions during adjustment (e.g. vibrations or air draught).
- Ensure that the scale has reached its operating temperature. ➔ *Chapter 18.1 'Technical data' on page 197*
- Only carry out the adjusting procedure with the standard weighing plate in place.
- During external adjusting, ensure that only the adjustment weight is on the weighing plate without any other load.
- Only handle adjustment weights with tweezers, weight handles or gloves. Avoid direct skin contact between your hands and adjustment weights.

6.3.1 Adjusting procedure with an external adjustment weight

To carry out the following steps, the scale must be in the setup menu mode.

1. ➤ Select **<AL>** → **<AL E H t>** (**[↑]** / **[↓]**) and confirm (**[→]**).

➔ The first applicable adjustment weight is displayed.

2. ➤ Select the required adjustment weight (**[↑]** / **[↓]**) and confirm (**[→]**).


(Information on adjustment weights ➔ *Chapter 18.1 'Technical data' on page 197*).


➔ **<Zero>** is displayed.

3. ➤ Make sure the weighing plate is free of any load and confirm by pressing **[→]**.

➔ **<P u t d>** → and the value of the adjustment weight are displayed.

4. ➤ Place the adjustment weight in the centre of the scale.

5. Confirm by pressing **[→]**.
 ➔ <BA It> → <rENULd> is displayed.


6. Remove the adjustment weight.
 ➔ The scale switches to operating mode.
 ➔ The adjusting process is complete.

If an adjustment error occurs (e.g. additional items on the weighing plate), the error message <ErronG> will appear. In this case, restart the scale and repeat the adjusting procedure.


6.3.2 Adjusting procedure with a user-defined adjustment weight

To carry out the following steps, the scale must be in the setup menu mode.

1. Select <CAL> → <CALEud> (**[↑]** / **[↓]**) and confirm (**[→]**).

➔ A numerical input field appears, and the active digit is displayed flashing.



2. Enter the weight value (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[→]**).

➔ <ZERO> is displayed.


3. Make sure the weighing plate is free of any load and confirm by pressing **[→]**.

➔ <PULld> → and the value of the adjustment weight are displayed.

4. Place the adjustment weight in the centre of the scale.

5. Confirm by pressing **[→]**.
 ➔ <BA It> → <rENULd> is displayed.

6. Remove the adjustment weight.

- ➡ The scale switches to operating mode.
- ➡ The adjusting process is complete.

If an adjustment error occurs (e.g. additional items on the weighing plate), the error message <ErOnG> will appear. In this case, restart the scale and repeat the adjusting procedure.

6.4 Setting the acoustic signal

<KEYtOn>

Settings

- <On>*: keypad sound enabled
- <OFF>*: keypad sound disabled

6.5 Setting the automatic switch-off function

<AutoFF>

Submenu	Description	Settings
<Mode>	Set switch-off mode	<ul style="list-style-type: none"> ■ <OFF>*: automatic switch-off function disabled ■ <Auto>: automatic switch-off after the time set at <Time>, provided there is no change of load ■ <only0>: automatic switch-off only when the display shows zero
<Time>	Set switch-off time	<ul style="list-style-type: none"> ■ <30s>* ■ <1m 00s> ■ <2m 00s> ■ <5m 00s> ■ <30m 00s> ■ <60m 00s>

6.6 Setting the automatic switch-on function

<AutoOn>

Settings

- <ON>: switch on the device automatically after connecting it to the mains voltage
- <OFF>*: disable the automatic switch-on function

6.7 Assigning of key functions

Certain scale keys can be programmed to perform different functions. This allows quick access to specific functions in operating mode.

<button>

Submenu	Description	Settings	
Key labelling (e.g. <CHANGE> ([CHANGE]) or <F-FKEY> ([F]), depending on the device)	Change the functionality of the selected key	<SPUSH>: short keystroke	<ul style="list-style-type: none"> ■ <toggle>: change the unit displayed ■ <unit>: set weighing unit ■ <showGt>: display gross and tare values ■ <hold>: execute hold function (only displayed when the "Weighing" application is active) ■ <mode>: change application ■ <Ptare>: open pre-tare settings ■ <chEcF>: When the "Counting" application is active: set the limits for piece counting When the "Tolerance weighing" application is active: open the settings menu for tolerances
		<LPUSH>: long keystroke	

Submenu	Description	Settings
		<ul style="list-style-type: none"> ■ <rEF>: set reference quantity (only displayed when the "Counting" application is active) ■ <OFF>: key disabled

6.8 Adjusting the backlight

<bL iGht>

Submenu	Description	Settings
<mode>	Set the lighting mode	<ul style="list-style-type: none"> ■ <ALWAYS>*: backlight always on ■ <t iNEr>: backlight switches off automatically after the time set at <t iNE>, provided there is no change of load ■ <noBL>: backlight off
<t iNE>	Set the lighting switch-off time	<ul style="list-style-type: none"> ■ <5>* ■ <30> ■ < 1n in> ■ <2n in> ■ <5n in> ■ <30n in>

6.9 Setting the taring range

<tArErG>

This allows the maximum taring range to be set between 10 % and 100 % (default value = 99 %). To do so, enter the taring range as a numerical value.

6.10 Setting the zero range (Zero Limit)

<ZL n t>

This allows the maximum zero range to be set between 1 % and 100 % (default value = 2 %). To do so, enter the zero range as a numerical value.

6.11 Setting automatic zero tracking

<Z E F A C F>

Settings

- <ON>*: automatic zero tracking ($\leq 3d$) enabled
- <OFF>: automatic zero tracking disabled

If small quantities of the weighing material are removed or added, the stability compensation may produce incorrect results (e.g. in cases of slow liquid outflow or evaporation). This function should therefore be deactivated for dosing operations involving minor weight fluctuations.

6.12 Setting the display speed / filter

<F I L T E R>

This allows the display speed to be set between 0 and 100 (0 = very fast for quiet environments, 100 = very slow for dynamic environments). To set the filter, enter a numerical value.

6.13 Setting the availability of display units

These settings allow you to determine which units are available in the application menu. For example, you can hide units that are not relevant to the application field.

<U N I T S>

Submenu	Description	Settings
For available units, see technical data	Selection of unit	<ul style="list-style-type: none"> ■ <ON>: make a unit available ■ <OFF>: lock unit

6.14 Setting the availability of applications

These settings allow you to determine which applications are available in the application menu. For example, you can hide applications that are not relevant to the application field.

Settings for an application can only be adjusted if it is not actively running in operating mode.

<ModE>

Submenu	Description	Settings
For available applications, see technical data	Selection of application	<ul style="list-style-type: none"> ■ <00>: make application available ■ <0FF>: lock application

6.15 Restoring factory settings

<RESE>

Once you have selected <RESE>, the message <SURE> will appear on the display. Once you have confirmed this message, the scale will switch off automatically. The next time you start the device, the settings will have been reset to their factory defaults. However, date and time will not be reset.

7 Interfaces and data output

7.1 Optional software

BalanceConnection

BalanceConnection enables communication between scales and peripheral devices. This software supports various interfaces and protocols to ensure flexible integration into existing systems. It is particularly useful for applications where centralised collection and processing of measurement data is required.

EasyTouch

With EasyTouch, your weighing data is sent directly to your PC, laptop or tablet for immediate further processing. EasyTouch also allows you to add other useful functions to your scale and tailor it perfectly to your specific field of application.

For further information about KERN software, please visit our website at ➔ www.kern-sohn.com.

7.2 Description of the interfaces

KUP (KERN Universal Port)

You can find all available KUP interface adapters in our online store at ➔ www.kern-sohn.com. Installation and setup are described in the relevant interface adapter documentation.

7.3 Setting up interfaces

7.3.1 Connecting wired peripheral devices

1. ➤ Switch off the scale and the peripheral device.
2. ➤ Disconnect the scale from the mains voltage.
3. ➤ Connect the KUP adapter to the scale.
4. ➤ Connect the peripheral device to the scale using a cable.
5. ➤ Reconnect the scale to the mains voltage.
6. ➤ Switch on the scale and the peripheral device.

➡ The interface is now ready for use. You may need to configure the interface parameters first. ➔ *Chapter 7.3.2 'Setting interface parameters' on page 150*

7.3.2 Setting interface parameters

<conf>

The settings options are only displayed once the interfaces have been fully initialised. During initialisation, the message <busy> is displayed. This may take a short while.

Only the settings for the connected interface are displayed.

Submenu	Submenu	Description	Settings
<rs232> or <usb-d>	<bAud>	Baud rate	<ul style="list-style-type: none"> ■ <600> ■ <1200> ■ <2400> ■ <4800> ■ <9600>* ■ <14400> ■ <19200> ■ <38400> ■ <57600> ■ <115200> ■ <230400>
	<dAtA>	Number of data bits	<ul style="list-style-type: none"> ■ <7db it> ■ <8db it>*
	<PAR it>	Parity	<ul style="list-style-type: none"> ■ <none>*: none ■ <odd>: odd ■ <EVEN>: even ■ <SPACE>: space ■ <PARF>: padding character
	<stop>	Number of stop bits	<ul style="list-style-type: none"> ■ <1st it>* ■ <2st it>
	<handsh>	Handshake	<ul style="list-style-type: none"> ■ <none>*

Submenu	Submenu	Description	Settings
	<Protoc>	Protocol	<FCP>*
<LAN>	<IP-Adr>	Display current IP address	
	<Subnet>	Display current subnet mask	
	<GtWAY>	Display current gateway address	
	<Reset>	Reset	

7.4 Configuring data output

7.4.1 Overview of output settings

<Print>

Submenu	Submenu	Sub-menu	Description	Settings
<IntFcE>			Select data output interface (only connected interfaces are shown)	<ul style="list-style-type: none"> ■ <r5232>* ■ <usb-d> ■ <ANa-oP> ■ <BLAN> ■ <LAN> ■ <BLE>
<Subn>			Enable/disable totalizing function	<ul style="list-style-type: none"> ■ <on>: enabled ■ <off>*: disabled

Submenu	Submenu	Sub-menu	Description	Settings
<PrNode>	<trIG>	<MANUAL>	Data output by pressing [PRINT]	<ul style="list-style-type: none"> ■ <on>*: enabled ■ <off>*: disabled
		<AutoPr>	Automatic data output if value is stable and positive	<div> <div><on>: enabled</div> <div><2rAnGE>: factor for d. This factor, when multiplied by d, gives the range of values that will not be printed.</div> </div>
				<off>*: disabled
		<cont>	Continuous data output (not available on verified models)	<div> <div><on>: enabled</div> <div> <div><intERU>: output interval in ms</div> <div><PZERO>: output when the scale is unloaded</div> <div><PSTab>: transfer stable values only (Cannot be changed on verified devices; set to <on> by default.)</div> </div> </div>
				<off>*: disabled
	<BEIGHt>	<SLPrt>	Output parameter: displayed value	<ul style="list-style-type: none"> ■ <on>: displayed weight value is output ■ <off>*: output disabled

Submenu	Submenu	Sub-menu	Description	Settings
		<Gross> <Net> <Tare>	Output parameters: Gross, Net, Tare (can be combined as required)	<Gross>: Gross ■ <on>*: enabled ■ <off>: disabled
				<Net>: Net ■ <on>*: enabled ■ <off>: disabled
				<Tare>: Tare ■ <on>*: enabled ■ <off>: disabled
				<Format>: format ■ <Long>: full designation ■ <Short>*: abbreviated designation
	<Layout>		Configure output layout	➔ <i>Chapter 7.4.3 'Print settings' on page 157</i>
	<Reset>		Reset output settings	■ <no>: do not reset settings ■ <yes>: reset settings

7.4.2 Setting output parameters

- Once an output mode has been enabled under **<Er IG>**, the other modes are automatically disabled. It is not possible to use multiple output modes at the same time.
- Configure the values to be output under **<PrNode>** → **<BE IGht>**.
→ Chapter 7.4.1 'Overview of output settings' on page 151

Manual data output **<MANUAL>**

If manual data output is selected, printing will only start after pressing **[PRINT]**.

1. Open the setup menu.
2. Select **<Pr int>** → **<PrNode>** → **<Er IG>** → **<MANUAL>** → **<on>** by pressing (**[↑]** / **[↓]**) and confirm your selection (**[→]**).
3. Press **[←]** repeatedly to exit the menu.
4. If necessary, place the weighing container on the scale and perform taring.
5. Place the load on the weighing plate.
6. Press **[PRINT]**.
➡ The result is displayed.

Automatic data output **<AutoPr>**

If automatic data output is enabled, printing starts as soon as the weight value is stable and positive.

In the settings for **<DrANGE>**, a factor for d can be defined to specify the range of values that will not be output.

Example

<DrANGE> = 5

d = 0.1 g

Result: values of 0.5 g or less will not be printed

1. Open the setup menu.

- 2. Select <Print> → <PrintMode> → <Unit> → <AutoPr> → <on> by pressing ([↑] / [↓]) and confirm your selection ([→]).
➡ <UnitAuto> is displayed.
- 3. Confirm by pressing [→].
➡ <onFF> is displayed.
- 4. Select the desired factor ([↑] / [↓]) and confirm ([→]).

<onFF>*	no factor
<1>	1 x d
<2>	2 x d
<3>	3 x d
<4>	4 x d
<5>	5 x d

- ➡ <UnitAuto> is displayed.
- 5. Press [←] repeatedly to exit the menu.
- 6. If necessary, place the weighing container on the scale and perform taring.
- 7. Place the load on the weighing plate.
➡ As soon as the "Stability" indicator (▲▲) appears, the result is displayed.

New output is only possible once the display has returned to zero and stabilised (remove the load from the scale, wait for the → zero display and for the "Stability" indicator (▲▲) to appear).

Continuous data output
<cont>

- 1. Open the setup menu.

2. Select <Pr int> → <PrModE> → <Er U> → <cont> → <on> by pressing ([↑] / [↓]) and confirm your selection ([→]).

➡ < intErU> is displayed.

3. Confirm by pressing [→].

➡ The display switches to numeric input.

➡ The most recent setting is displayed flashing

4. Enter the interval in milliseconds (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

➡ The display switches to < intErU>.

5. Select <PZEro> by pressing ([↑] / [↓]) and confirm ([→]).

6. Select the desired setting ([↑] / [↓]) and confirm ([→]).

<on>	continuous output of 0 values (unloaded scale)
<off>	no output of 0 values (unloaded scale)

➡ The display switches to <PZEro>.

7. Press [←] repeatedly to exit the menu.

8. If necessary, place the weighing container on the scale and perform taring.

9. Place the load on the weighing plate.

➡ The weight values are output at defined intervals.

7.4.3 Print settings

<Print> → <Layout>

Submenu		Description	Settings
<USER>	<nonE>	Enable / disable default layout	<ul style="list-style-type: none"> ■ <ON>*: enable default layout (other layouts are automatically disabled) ■ <OFF>: disable default layout (default layout is automatically disabled when another layout is enabled)
	<MODEL>	Print the model name of the scale	<ul style="list-style-type: none"> ■ <OFF>*: disable custom layout setting (custom layout is automatically disabled when another layout is enabled) ■ <ON>: enable custom layout setting (other layouts are automatically disabled)
	<SERIAL>	Print the scale's serial number	
	<AL ID>	Print alibi ID (only for scales with an alibi memory)	
	<DATE>	Print the date (only for scales with an alibi memory)	
	<TIME>	Print the time (only for scales with an alibi memory)	

Print layouts (examples)

The print output depends on the following settings:

- Settings in <UNIT>, <WEIGHT> and <Layout>
- Selected display unit

The wide range of setting options allows for a variety of possible print layouts. The following sections provide some examples of settings and print layouts.

Print layout: standard

Settings ■ Standard unit ■ <GntPrt> - <Gross> → <on> - <Net> → <on> - <Tare> → <on> - <Format> → <Short>			
N: 500.0 g T: 0.0 g G: 500.0 g <i>Fig. 2: Simple weighing</i>	N: 500.0 g T: 0.0 g G: 500.0 g PER: 100.0 % <i>Fig. 3: Percentage weighing</i>	N: 500.0 g T: 0.0 g G: 500.0 g PCS: 5 pcs UW: 100.0 g REF: 5 pcs <i>Fig. 4: Counting</i>	

Settings ■ Standard unit ■ <GntPrt> - <Gross> → <on> - <Net> → <on> - <Tare> → <on> - <Format> → <Long>			
Net weight: 500.0 g Tare weight: 0.0 g Gross weight: 500.0 g <i>Fig. 5: Simple weighing</i>	Net weight: 500.0 g Tare weight: 0.0 g Gross weight: 500.0 g FFA: 550.0 ffa <i>Fig. 6: Multiplication (FFA)</i>		

Print layout: single value (net value)

Settings

- Standard unit
- <SGLPrt> → <on>

	500.0 g
PCS:	5 pcs
UW:	100.0 g
REF:	5 pcs
	500.0 g
PER:	100.0 %

Fig. 7: Single value (designation varies depending on the unit setting)

Print layout: custom

Model No: PCJ 6000-1M
 Serial No: N/A
500.1 g

Fig. 8: Custom printout

The model and serial number can be selected or deselected individually. "N/A" appears if no serial number has been set.

The **text in bold** can be customized by changing the settings for <BE iGhE> and unit.

Print layout: totalization

Layout	<uBEr>	<nonE>
Header	Model No: PCJ 6000-1M Serial No: N/A	-
Data	No. = 1 Net weight: 0.0 g Gross weight: 0.0 g TOTAL = 0.0 g No. = 2 Net weight: 0.0 g Gross weight: 0.0 g TOTAL = 0.0 g	
Total	No. = 3 TOTAL = 0.0 g	
Footer	-	

7.5 KCP - KERN Communications Protocol (KERN interface protocol)

KCP is a standardised interface command set for KERN scales, which allows users to retrieve and control a wide range of parameters and device functions. This makes it very easy to connect KERN devices with KCP to computers, industrial control systems and other digital systems. A detailed description can be found in the “KERN Communications Protocol” manual, available in the download section of our KERN website (→ www.kern-sohn.com).

KCP is based on simple ASCII commands and responses. Each interaction consists of a command, which may include arguments separated by spaces, and ends with <CR> <LF>.

You can query the KCP commands supported by your scale by sending the command „I0“ followed by <CR><LF>.

Tab. 1: Selection of commonly used KCP commands

Command	Description
I0	Show all implemented KCP commands
S	Send a stable value
SI	Send current value (even if unstable)
SIR	Send the current value (even if unstable) and repeat
T	Tare
Z	Zero setting

Example: possible responses to the command "S"

Possible answer	Description
S_S_100.00_g	The command has been executed successfully. A stable weight value is returned.
S_I	The menu is open, another command is currently being executed, or a timeout has been reached.
S_+ or S_-	Overload or underload

8 Basic operation

8.1 Switching on the device

WARNING

Serious injuries due to damaged power supply units

Damaged power supply units can lead to short circuits, fires and serious injuries.

- Before each use, check the power supply unit and mains cable for visible damage. Do not use the power supply unit if it is damaged.

WARNING

Injuries due to incorrect handling of power supply units

Incorrect handling of power supply units can lead to short circuits, fires and serious injuries.


- Observe the safety precautions for handling power supply units. ➔ *'Safe handling of power supply units and power supplies' on page 117*
- Avoid moisture at the mains connection, power supply unit and socket.
- Only connect the device to a properly installed socket.
- Make sure that the power supply unit cable is never kinked or pinched.
- Ensure that the power supply unit is accessible at all times.
- Ensure that the mains cable does not pose a tripping hazard.
- Only connect the power supply unit to the mains if the information on the device (sticker) and the mains voltage at the place of use are identical.

NOTICE

Incorrect weighing results

Ensure that the scale has reached its operating temperature before you start weighing. ➔ [Chapter 18.1 'Technical data' on page 197](#)

- ▶ Press **[ON-OFF]**.


- ➔ The display lights up.
- ➔ The scale is ready for operation.

8.2 Switching off the device

NOTICE

Do not leave a load on the weighing plate if you are planning to switch off the scale for a longer period of time. This could damage the measuring element.

- ▶ Press and hold **[ON-OFF]**.


- ➔ The display switches off.
- ➔ The scale is switched off.

8.3 Zero setting


Zero setting a scale sets the display to zero so that no weight is displayed when the weighing plate is unloaded. This ensures measurement accuracy by compensating for slight deviations or external influences.

Zero setting is only possible within +/- 2 % of the maximum weighing range. An error message appears if the value is above or below this range limit. ➔ [Chapter 17.1 'Error messages' on page 194](#)

1. ➤ Unload the scale.

2. ➤ Press **[TARE]**.

- ➔ Zero and the "Zero display" indicator (>0<) are displayed.
- ➔ Scale is set to zero.

8.4 Taring


At the touch of a button, the tare weight of any weighing container can be compensated so only the net weight of the weighed material will be displayed in subsequent weighing operations.

- After removing the container, the weight of the container appears as a negative value.
- To delete the stored tare value, remove the load from the weighing plate and tare or zero the unloaded scale.
- The taring process can be repeated any number of times, for example when weighing several components to form a mixture (add-on weighing). The limit is reached when the taring range is fully utilised.
- With the PRE-TARE function, a tare weight can be saved permanently. ➔ [Chapter 8.4.2 'Saving a known tare weight \(pre-tare\)' on page 164](#)


8.4.1 Tare

Tare


1. ➤ Place the empty container on the weighing plate.

- ➔ The weight of the container and the "Stability" indicator (▴▾) are displayed.

Weighing with tare

2. Press **[TARE]**.

- ➡ The zero display and the "Net weight" indicator (**NET**) are displayed.
- ➡ The weight of the container is stored as a tare.

3. Fill the material to be weighed into the container.

- ➡ The weight value and the "Stability" indicator (▴▴) are displayed.
- ➡ The result can be read.



8.4.2 Saving a known tare weight (pre-tare)

The pre-tare function allows the weight of a container to be permanently stored in the scale. This is helpful in processes where containers of the same weight are used.

- The tare weight remains valid until a new tare weight is saved or the tare weight is deleted.
- A previously stored tare weight can be recalled: open the application menu and select <P**RE**TARE> → <R**EC**ALL>.
- There are two ways to delete the tare weight:
 - Unload the scale and press **[TARE]**.
 - Open the application menu and go to <P**RE**TARE> → <C**LEAR**>.


Set and save tare weight

1. Place the empty container on the weighing plate.

2. Press and hold **[TARE]**.


- ➡ The application menu opens.

3. Select <PRETARE> → <ACTUAL> by pressing [↑] / [↓] and confirm ([→]).


- ➔ The zero display and the "Pre-tare" indicator (**PRETARE**) as well as the "Net weight" indicator (**NET**) are displayed.
- ➔ The tare weight has been saved.

Enter tare weight manually


1. Press and hold [TARE].

- ➔ The application menu opens.

2. Select <PRETARE> → <MANUAL> by pressing [↑] / [↓] and confirm ([→]).

- ➔ The display switches to numeric input.


3. Enter the tare weight (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

- ➔ The "Pre-tare" indicator (**PRETARE**) and the "Net weight" indicator (**NET**) are displayed.
- ➔ The tare weight has been saved and will be displayed as a negative value.
- ➔ Weighing can now be carried out using the appropriate container.

8.5 Selecting an application


Open the application menu

1. Press and hold [TARE].

- ➔ <APcNEr> followed by the currently active application is displayed.
- ➔ The scale is in application menu mode.

Select application

2. Select <NODE> ([↑] / [↓]) and confirm ([→]).

- ➔ The available applications are displayed.

3. → Select the desired application (**[↑]** / **[↓]**) and confirm (**[→]**).

<WE ih>	Weighing <i>➔ Chapter 9 'Weighing application' on page 167</i>
<count>	Counting <i>➔ Chapter 10 'Counting application' on page 172</i>
<chEcH>	Tolerance weighing <i>➔ Chapter 11 'Tolerance weighing application' on page 181</i>

- ➔ The scale returns to operating mode.
- ➔ The application is now active.

Once an application has been selected, the device returns to operating mode. You must then open the application menu again to view the application-specific settings.

9 Weighing application


Instructions for zero setting and taring are provided in [Chapter 8 'Basic operation' on page 161](#). Further specific setting options can be found in the following chapters.

9.1 Overview of the "Weighing" application menu

Level 1	Description	Further information
<P <small>RE</small> T <small>ARE</small> >	Save known tare weight (pre-tare)	Chapter 8.4.2 'Saving a known tare weight (pre-tare)' on page 164
<h <small>OLD</small> >	Hold reading (data hold)	Chapter 9.4 'Hold function' on page 168
<u <small>N</small> I <small>T</small> >	Set weighing units	Chapter 6.13 'Setting the availability of display units' on page 147
<M <small>ODE</small> >	Select application	Chapter 8.5 'Selecting an application' on page 165

9.2 Simple weighing

To use this function, the "Weighing" application must be active.

1. Check whether the display reads zero and, if necessary, reset it to zero using **[TARE]**.
2. Place the load on the weighing plate.
 - ➔ The weight value and the "Stability" indicator (▲▲) are displayed.
 - ➔ The result can be read.

9.3 Underfloor weighing

Underfloor weighing allows objects that cannot be placed on the weighing plate due to their size or shape to be weighed.

CAUTION

Risk of breakage due to overloading of the hook

Falling loads may cause injury.

- Never exceed the scale's specified maximum load (Max.).
- Make sure there are no living creatures or objects under the load that could be damaged.

NOTICE

Damage to the scale caused by dirt getting inside

Dirt (e.g. dust) entering the measuring system can cause damage and lead to incorrect results.

- Close the opening for the weigh-below hook once underfloor weighing is complete (to prevent dust ingress).

1. ➤ Switch off the scale.
2. ➤ Turn the scale over.
3. ➤ Open the cover on the bottom side of the scale.
➔ *Chapter 3.2 'Overview of the device' on page 123*
4. ➤ Place the scale over an opening.
5. ➤ Screw the hook in tightly.
6. ➤ If necessary, attach the weighing container and perform taring.
7. ➤ Attach the load and weigh it.

9.4 Hold function

The hold function allows you to 'freeze' the displayed weight. Once the load is removed from the scale, the weight value remains visible on the display for 10 seconds.

To use this function, the "Weighing" application must be active.

-
1. ➤ Press and hold **[TARE]**.
➡ The application menu opens.
 2. ➤ Select **<hold>** by pressing **[↑]** / **[↓]**, but do not confirm yet.
 3. ➤ Place the load on the scale.
 4. ➤ Confirm by pressing **[→]**.
➡ The first stable weight reading is 'frozen' and indicated by the "Hold" symbol (**HOLD**).




For quicker access, the hold function can be assigned to a key.

9.5 Percentage weighing

The "Percentage weighing" function allows you to calculate the percentage of the sample weight in relation to a reference weight.

To use this function, the "Weighing" application must be active.

-
1. ➤ Check whether the display reads zero and, if necessary, reset it to zero using **[TARE]**.
 2. ➤ Press and hold **[TARE]**.
➡ The application menu opens.

3. Select $\langle \text{unit} \rangle$ by pressing $\uparrow$ / $\downarrow$ and confirm ($\rightarrow$).

➔ The display switches to the unit selection menu.

4. Place the reference weight, which corresponds to 100 %, on the scale.

5. Select "%" by pressing $\uparrow$ / $\downarrow$ and confirm ($\rightarrow$).

➔ The weight value of the reference weight will be displayed flashing.

6. Press $\rightarrow$.




➔ The scale switches to operating mode and displays 100 %.

➔ When weighing additional loads, their weight is displayed as a percentage of the reference weight.

9.6 Weighing with a multiplication factor

This function allows the weight value to be multiplied by any factor. This means, for example, that a known error factor can be taken into account immediately when determining the weight.

To use this function, the "Weighing" application must be active.

1. Press and hold **TARE**.


➔ The application menu opens.

2. Select $\langle \text{unit} \rangle$ by pressing $\uparrow$ / $\downarrow$ and confirm ($\rightarrow$).

➔ The display switches to the unit selection menu.

3. Select "FFA" by pressing $\uparrow$ / $\downarrow$ and confirm ($\rightarrow$).

➔ The display switches to numeric input.

4. Enter the factor (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[→]**).

➡ The scale switches to operating mode and displays the weight multiplied by the factor.

10 Counting application

The Counting application allows you to count several parts placed on the weighing plate.

10.1 Overview of the "Counting" application menu

Level 1	Level 2	Level 3	Description	Further information
<REF>			Set the reference quantity	➔ Chapter 10.2 'Setting the reference quantity' on page 174
<PRE>			Save known tare weight (pre-tare)	➔ Chapter 8.4.2 'Saving a known tare weight (pre-tare)' on page 164
<CHECK>	<TYPE>	<AB>	Check quantity (limit values)	➔ Chapter 10.5 'Quantity checking (limit values)' on page 179
		<REL>	Weigh the target number of items (deviations)	➔ Chapter 10.4 'Weighing the target number of items' on page 177
	With setting <TYPE> = <AB>			
	<L INPP>		Upper limit	
	<L INLB>		Lower limit	
	With setting <TYPE> = <REL>			
	<VALUE>		Target number of items	

Level 1	Level 2	Level 3	Description	Further information
	<ErruPP>		Upper tolerance	
	<ErrLoB>		Lower tolerance	
	<noB iG>		Define a range to be excluded from the tolerance assessment	
	<bEEPER>	<ch-off> (value within the tolerances / limits)	Sound when tolerance is exceeded ■ <off>* : no sound	
		<ch-Lo> (below the lower tolerance / below the lower limit)	■ <LoB> : slow beeps ■ <LoD>: normal beeps	
		<ch-h i> (upper tolerance / upper limit exceeded)	■ <FASt> : rapid beeps ■ <cont> : continuous beep	
	<rESEt>		Reset values	
<un it>			Set weighing units	➔ Chapter 6.13 'Setting the availability of display units' on page 147


Level 1	Level 2	Level 3	Description	Further information
<ModE>			Select application	➔ Chapter 8.5 'Selecting an application' on page 165

10.2 Setting the reference quantity


Before the scale can determine the number of items, it needs to calculate the average weight per item.

- The higher the reference quantity, the greater the counting accuracy.
- For small items or items that vary significantly one from another, the reference value must be set particularly high.
- The minimum part weight for piece counting is specified in the technical data.

Option 1: Select a predefined reference quantity

1. (1) If necessary, place the weighing container on the scale and perform taring.


2. (2) Load the desired number of reference items.

3. (3) Confirm the loaded reference quantity by pressing the corresponding button (5x, 10x, 20x).

➔ The scale calculates the average unit weight and then displays the number of items.

➔ (4) The scale is now ready for piece counting.

Option 2: Set a custom reference quantity

1. (1) If necessary, place the weighing container on the scale and perform taring.

2. (2) Load the desired number of reference items.

3. (3) Press and hold **[10x/Ref n]**.


➔ The display switches to numeric input.

➔ The most recent setting is displayed flashing





4. Enter the required reference quantity (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[↵]**).

➔ The scale returns to operating mode and displays the number of items.

➔ (4) The scale is now ready for piece counting.




Option 3: Enter the average piece weight

1. Press and hold **[TARE]**.
 - ➔ The application menu opens.
2. Select **< rEF > → < inPut >** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
 - ➔ The display switches to the weighing unit settings.
3. Select the weighing unit (**[↑]** / **[↓]**) and confirm (**[→]**).
 - ➔ The display switches to the decimal place settings.
4. Select the decimal place (**[↑]** / **[↓]**) and confirm (**[→]**).
 - ➔ The display switches to numeric input.
5. Enter the piece weight (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[→]**).
 - ➔ The scale returns to operating mode and displays the number of items.
 - ➔ The scale is ready for piece counting.

10.3 Simple piece counting


1. Set the reference quantity. ➔ [Chapter 10.2 'Setting the reference quantity' on page 174](#)
2. If necessary, place the weighing container on the scale and perform taring.
3. Weigh the pieces.
 - ➔ The quantity and the "Stability" indicator (▲▲) are displayed
 - ➔ The number of items can be read.

Pressing **[CHANGE]** toggles the display between the current quantity and the total weight (default setting).


10.4 Weighing the target number of items

Select application mode


1. Press and hold **[TARE]**.
➔ The application menu opens.
2. Select **<chEcK>** → **<TYPE>** → **<REL>** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
➔ Application mode is active.

Set target quantity and tolerances

1. Set the reference quantity. ➔ [Chapter 10.2 'Setting the reference quantity' on page 174](#)
2. Press and hold **[TARE]**.
➔ The application menu opens.

Enter target quantity


3. Select **<chEcK>** → **<VALUE>** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
➔ The display switches to numeric input.


4. Enter the target quantity (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[→]**).
➔ The display switches to **<VALUE>**.

Enter upper tolerance

5. Select **<ErruPP>** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
➔ The display switches to the unit selection menu.

Counting application


6. Select the desired unit ([↑] / [↓]) and confirm ([→]).
 ➔ The display switches to numeric input.

7. Enter the upper tolerance (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).
 ➔ The display switches to <ErruPP>.

Enter lower tolerance


8. Select <ErrLoB> by pressing [↑] / [↓] and confirm ([→]).
 ➔ The display switches to the unit selection menu.

9. Select the desired unit ([↑] / [↓]) and confirm ([→]).
 ➔ The display switches to numeric input.



10. Enter the lower tolerance (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).
 ➔ The display switches to <ErrLoB>.

11. Press [←] repeatedly to exit the menu.

Weigh target number of items

1. If necessary, place the weighing container on the scale and perform taring.

2. Weigh the pieces.
 ➔ Quantity, "Stability" indicator (▲▲) and tolerance marks are displayed.
 ➔ The result can be read.

Tolerance mark	Description
	Target quantity exceeded
	Number of items within the tolerances
	Target quantity not reached

- Pressing **[CHANGE]** toggles the display between the current quantity and the total weight (default setting).
- In addition to the tolerance marks, an audible signal can also be set to check the limits. Further information on this can be found in the application menu.
- To delete the values, proceed as follows: Open the application menu and go to **<CHECK> → <RESET>**.

10.5 Quantity checking (limit values)

APCNEr

1. Press and hold **[TARE]**.
➔ The application menu opens.
2. Select **<CHECK> → <TYPE> → <ABS>** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
➔ Application mode is active.
3. Set the reference quantity. ➔ *Chapter 10.2 'Setting the reference quantity' on page 174*
4. Select **<CHECK> → <L. UPPER>** by pressing **[↑]** / **[↓]** and confirm (**[→]**).
➔ The display switches to numeric input.
5. Enter the upper limit (select a number: **[↑]** / **[↓]**, next number: **[→]**) and confirm (**[→]**).
➔ The display switches to **<L. UPPER>**.

Enter the upper limit

L. UPPER


000000

000079

L. UPPER

Enter the lower limit

Check the number of items

6. Select <L nLoB> by pressing [↑] / [↓] and confirm ([→]).

➡ The display switches to numeric input.

7. Enter the lower limit (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

➡ The display switches to <L nLoB>.

8. Press [←] repeatedly to exit the menu.

1. If necessary, place the weighing container on the scale and perform taring.

2. Weigh the pieces.

➡ Quantity, "Stability" indicator (▲▲) and tolerance marks are displayed.

➡ The result can be read.

Tolerance mark	Description
↑ HI	Upper limit exceeded
OK	Weight within the limits
↓ LO	Below the lower limit

- In addition to the tolerance marks, an audible signal can also be set to check the limits. Further information on this can be found in the application menu.
- To delete the values, proceed as follows:
Open the application menu and go to <chEcF> → <rE5Et>.

11 Tolerance weighing application

The tolerance weighing application allows you to quickly check whether a weight value falls within certain limits.

11.1 Overview of the "Tolerance weighing" application menu

Level 1	Level 2	Level 3	Description	Further information
<chEcH>	<TYPE>	<AB5>	Check weight (upper and lower limits)	➔ Chapter 11.3 'Weight checking (limit values)' on page 184
		<rEL>	Weigh target weight (deviations)	➔ Chapter 11.2 'Weighing the target weight' on page 182
	With setting <TYPE> = <AB5>			
	<L nUPP>		Upper limit	
	<L nLOB>		Lower limit	
	With setting <TYPE> = <rEL>			
	<UALUE>		Target weight	
	<ErrUPP>		Upper tolerance	
	<ErrLOB>		Lower tolerance	
	<nob iG>		Define a range to be excluded from the tolerance assessment	




Level 1	Level 2	Level 3	Description	Further information
	<bEEPER>	<ch-off> (value within the tolerances / limits)	Sound when tolerance is exceeded	
		<ch-Lo> (below the lower tolerance / below the lower limit)	<ul style="list-style-type: none"> ■ <off>* : no sound ■ <slow> : slow beeps 	
		<ch-Hi> (upper tolerance / upper limit exceeded)	<ul style="list-style-type: none"> ■ <std>: normal beeps ■ <fast>: rapid beeps ■ <cont>: continuous beep 	
	<reset>		Reset values	
<PtarE>			Save known tare weight (pre-tare)	➔ Chapter 8.4.2 'Saving a known tare weight (pre-tare)' on page 164
<unit>			Set weighing units	➔ Chapter 6.13 'Setting the availability of display units' on page 147
<mode>			Select application	➔ Chapter 8.5 'Selecting an application' on page 165

11.2 Weighing the target weight

APcNEr

- ➔ Press and hold **[TARE]**.
 - ➔ The application menu opens.
- ➔ Select <chEcH> → <type> → <rel> by pressing **[↑]** / **[↓]** and confirm (**[→]**).
 - ➔ Application mode is active.

Enter target weight**Enter upper tolerance****Enter lower tolerance****Weigh target weight**

3. Select <CHECK> → <VALUE> by pressing [↑] / [↓] and confirm ([→]).

➔ The display switches to numeric input.

4. Enter your target weight (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

➔ The display switches to <VALUE>.

5. Select <ErruPP> by pressing [↑] / [↓] and confirm ([→]).

➔ The display switches to the unit selection menu.

6. Select the desired unit ([↑] / [↓]) and confirm ([→]).

➔ The display switches to numeric input.

7. Enter the upper tolerance (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

➔ The display switches to <ErruPP>.

8. Select <ErrLoB> by pressing [↑] / [↓] and confirm ([→]).

➔ The display switches to the unit selection menu.

9. Select the desired unit ([↑] / [↓]) and confirm ([→]).

➔ The display switches to numeric input.

10. Enter the lower tolerance (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

➔ The display switches to <ErrLoB>.

11. Press [←] repeatedly to exit the menu.

1. If necessary, place the weighing container on the scale and perform taring.

2. Weigh the material to be weighed.

- ➔ Weight, "Stability" indicator (▲▲) and tolerance marks are displayed.
- ➔ The result can be read.

Tolerance mark	Description
	Upper tolerance limit exceeded
	Weight within tolerance limits
	Below the lower tolerance limit

- In addition to the tolerance marks, an audible signal can also be set to check the limits. Further information on this can be found in the application menu.
- To delete the values, proceed as follows:
Open the application menu and go to <CHECK> → <RESET>.

11.3 Weight checking (limit values)

1. Press and hold [TARE].

- ➔ The application menu opens.


2. Select <CHECK> → <TYPE> → <Ab> by pressing [↑] / [↓] and confirm ([→]).

- ➔ Application mode is active.

3. Select <CHECK> → <L 100PP> by pressing [↑] / [↓] and confirm ([→]).

- ➔ The display switches to numeric input.


Enter the upper limit

4. Enter the upper limit (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).

- ➔ The display switches to <L 100PP>.

Enter the lower limit

Check weight

5. Select <L nLoB> by pressing [↑] / [↓] and confirm ([→]).
➔ The display switches to numeric input.
6. Enter the lower limit (select a number: [↑] / [↓], next number: [→]) and confirm ([→]).
➔ The display switches to <L nLoB>.
7. Press [←] repeatedly to exit the menu.
1. If necessary, place the weighing container on the scale and perform taring.
2. Weigh the material.
➔ Weight, "Stability" indicator (▴▴) and tolerance marks are displayed.
➔ The result can be read.

Tolerance mark	Description
HI	Upper limit exceeded
OK	Weight within the limits
LO	Below the lower limit

- In addition to the tolerance marks, an audible signal can also be set to check the limits. Further information on this can be found in the application menu.
- To delete the values, proceed as follows: Open the application menu and go to <chEcK> → <rESEt>.

12 Other features

12.1 Applying the totalizing function

The totalizing function allows you to add together the weight values from individual weighings.

To use the totalizing function, your scale must be connected to a printer or a computer.

Enable the function

1. Open the setup menu.
2. Select $\langle \text{Print} \rangle \rightarrow \langle \text{Sum} \rangle \rightarrow \langle \text{On} \rangle$ by pressing $\uparrow$ / $\downarrow$ and confirm ($\rightarrow$).
3. Press $\leftarrow$ repeatedly to exit the menu.

Apply the totalizing function

1. If necessary, place the weighing container on the scale and perform taring.
2. Place the first load to be weighed on the weighing plate.
 - The "Stability" indicator (▴ ▾) appears.
3. Press **[PRINT]**.
 - The value is stored in the total memory and the display switches to $\langle \text{Sum} \rangle$, followed by the current weight value.
 - The "Total" indicator (Σ) appears.
 - The result is sent to a printer or PC.
4. Remove the weighed material.
5. Add any further loads to be weighed, as in steps 2 to 4.

Display total weight and clear the total memory

- The number of weighings is displayed via $\langle \text{Sum} \rangle$, $\langle \text{Sum}2 \rangle$, $\langle \text{Sum}3 \rangle$ and so on.
- ▶ Press and hold **[PRINT]**.
 - The number of weighings and the total weight are displayed.
 - The total memory is cleared and the "Total" indicator (Σ) disappears.
 - The result is sent to a printer or PC.

13 Care, maintenance and servicing

13.1 Cleaning

13.1.1 Safety during cleaning

DANGER

Electric shock due to live parts

Contact with live parts can lead to electric shock and burns.

- Installation, cleaning, maintenance and servicing work may only be carried out by qualified personnel.
- Disconnect live parts from the power supply.
- Only carry out installation, cleaning, servicing and maintenance work on devices that are switched off and disconnected from the power supply, unless otherwise specified in the relevant instructions.
- Follow the instructions in this document.

NOTICE

Material damage due to incorrect cleaning

Unsuitable cleaning agents can cause damage to the device.

- Do not use any aggressive cleaning agents (e.g. solvents).
- Do not use cleaning agents containing caustic soda, acetic, hydrochloric, sulphuric or citric acid on stainless steel parts.
- Do not use metal brushes or cleaning sponges made of steel wool.
- Do not open the device. There are no parts inside that may be cleaned.

13.1.2 Recommended cleaning

- 1.** ➤ Switch off the device.
- 2.** ➤ Disconnect the device from the power supply.

3. ➤ Immediately remove any spilled material (e.g. loose sample residues or powder) with a brush or a handheld vacuum cleaner.
4. ➤ Clean the device with a gentle cleaning agent (e.g. soap suds) and a soft cloth.
5. ➤ Use a dry, lint-free cloth to wipe the device.

13.2 Servicing

13.2.1 Safety during servicing

DANGER

Electric shock due to live parts

Contact with live parts can lead to electric shock and burns.

- Installation, cleaning, maintenance and servicing work may only be carried out by qualified personnel.
- Disconnect live parts from the power supply.
- Only carry out installation, cleaning, servicing and maintenance work on devices that are switched off and disconnected from the power supply, unless otherwise specified in the relevant instructions.
- Follow the instructions in this document.

13.2.2 Maintenance schedule

Interval	Servicing activity	Description
Regularly	Levelling ➔ <i>Chapter 5.2.4 'Levelling the device' on page 133</i>	<ul style="list-style-type: none">■ Level the scale after moving it to a new location.■ Check the levelling regularly.
If necessary	Charge battery	Charge the rechargeable battery when it is flat.

Interval	Servicing activity	Description
Regularly / after relocating the device	Adjusting ➔ <i>Chapter 6.3 'Adjustment settings' on page 141</i>	Adjust the scale at regular intervals and after each change of location. KERN & SOHN offers test weights for this purpose.

13.3 Control of inspection, measuring and test equipment

As part of quality assurance, the metrological characteristics of the scale and any test weights must be checked at regular intervals. The responsible user must define a suitable interval for this, as well as the nature and scope of the test. Information regarding the monitoring of scales and the test weights required for this purpose is available on the KERN website (➔ www.kern-sohn.com). KERN's accredited Calibration Laboratory offers a fast and cost-effective service for the calibration of test weights and scales (traceability to the national measurement standard).

14 Decommissioning and storage

14.1 Taking the device out of operation

1. ➤ When operating with a rechargeable battery: To prevent deep discharge, charge the battery if the device is to be stored for a longer period of time.
2. ➤ Switch off the device.
3. ➤ Disconnect the device from the power supply:
 - Disconnect the power supply from the mains and unplug it from the device.
 - Remove the rechargeable battery from the device.
4. ➤ Remove all loose / moving parts from the device.
5. ➤ Remove all interface cables.
6. ➤ Pack the device and all accompanying parts neatly in the original packaging.
 - ➡ The device can now be stored ➔ [Chapter 14.2 'Storing the device' on page 190](#).

14.2 Storing the device

NOTICE

Damage caused by incorrect storage

Improper storage can damage the device.

- Take the device out of operation according to the instructions in this document ➔ [Chapter 14.1 'Taking the device out of operation' on page 190](#).
- Store the device in the original packaging.
- Do not store the device outdoors.
- Observe the environmental conditions in the technical data.

15 Warranty and goods return

15.1 Warranty

The warranty claim only applies to defects that were already present at the time of purchase. Damages or defects caused by improper use, incorrect handling or normal wear and tear are excluded from the warranty.

15.2 Goods return

- *A return is not possible if you have used the device improperly and it has been damaged as a result.*
- *The original packaging, including all associated parts, must be used for returns.*
- *In cases of returns, the complete scope of delivery must be sent back.*

Procedure for returns

- 1.** ➤ Check that all parts of the delivery are present.
- 2.** ➤ Remove all loose / moving parts and connected cables from the device.
- 3.** ➤ Reattach the transport safety lock.
- 4.** ➤ Reapply any protective film that may be present.
- 5.** ➤ Pack all parts of the delivery in the original packaging and secure them so that they cannot slip and are protected from damage.

16 Disposal

ENVIRONMENT

Damage to people and the environment due to improper disposal

Non-degradable materials (e.g. plastics or metals) can be released into the environment and cause damage if disposed of improperly.

Materials that are not recycled are a waste of resources and lead to increased environmental pollution.

Electrical devices may contain substances that are not biodegradable and are harmful to the environment and human health.

Batteries and rechargeable batteries may contain substances that are harmful to the environment and human health.

- Dispose of the device via appropriate recycling measures or recycling centres in accordance with the applicable national or regional law of the place of use.

Disposal of electronic devices

The crossed-out wheeled bin symbol indicates that waste from electronic equipment must not be disposed of with household waste. Disposal must be arranged by the operator in accordance with the applicable national or regional law of the place of use.

If this device contains confidential data, it is the responsibility of the operator to delete this data before disposal.

Remove batteries and rechargeable batteries before disposal of the device and dispose of them separately ➔ *'Disposal of batteries and rechargeable batteries' on page 192.*

Disposal of batteries and rechargeable batteries

The crossed-out wheeled bin symbol indicates that batteries and rechargeable batteries must not be disposed of with household waste. Disposal of batteries and rechargeable batteries must be arranged by the operator in accordance with the applicable national or regional law of the place of use.

We are obliged to inform end users in accordance with the EU Battery regulation and the applicable national regulations:

- End users are legally obliged to return all used (rechargeable) batteries.
- Batteries and rechargeable batteries can be returned free of charge after use to municipal collection centres or retailers. Non-used batteries/rechargeable batteries must be protected against short circuits.
- Returns to us are limited to usual household quantities of batteries and rechargeable batteries that we supply or have supplied.
- If batteries or rechargeable batteries contain cadmium (Cd), mercury (Hg) or lead (Pb), the respective symbol is shown below the waste bin symbol.

17 Troubleshooting

17.1 Error messages

Fault description	Cause	Remedy
0L 00 00	Zero-setting range exceeded	Ensure that the zero-setting range is within +/- 2% of the maximum weight
undEr0	Value below the zero-setting range	
instAb	Load unstable	<ul style="list-style-type: none"> ■ Place the load gently ■ Create stable environmental conditions ■ Eliminate electromagnetic fields
Er000	Adjusting error	Perform the adjusting again
5Etrtc	Date and time are incorrect or have not been set	Set date and time in the setup menu
0unAU	Alibi memory not available or faulty	Contact the KERN technical service department
no232	RS232 not available	Connect the RS232 KUP
noBLAn	Wi-Fi not available	Connect the Wi-Fi KUP
L---	Underload	The weight must be at least equal to the minimum weight
----	Overload	Apply less weight
LoBAE	Battery low	Change battery / Charge rechargeable battery

17.2 Troubleshooting table

Fault description	Cause	Remedy
Display does not switch on	Power supply not connected	Connect the device to the power supply
	Device is not switched on	Switch on the device
	Display defective	Contact the KERN technical service department

Fault description	Cause	Remedy
Device cannot be switched on	Power key defective	Contact the KERN technical service department
	Power supply unit defective or not connected	Check and plug in the power supply unit
	Mains power outage	Wait
	Rechargeable battery empty	Charge battery
	Rechargeable battery not connected	Connect rechargeable battery
Display does not respond to load changes	Display defective	Contact the KERN technical service department
	Load cell defective	Contact the KERN technical service department
	Transport safety lock not removed	Remove the transport safety lock ➔ Further information on page 132
Displayed weight is constantly changing	Load is moving	Place the load gently
	Unstable environmental conditions (e.g. draught / air movement, vibrations of the table or floor)	Stabilise the environmental conditions or move the device to another location
	Weighing plate is in contact with foreign objects	Remove foreign objects
	Electromagnetic fields / Electrostatic charge	If possible, switch off the interfering device or change the location of the scale
	Load cell defective	Contact the KERN technical service department
Result is incorrect	Scale is off zero	Unload the scale and set it to zero
	Adjustment is no longer correct	Carry out adjusting
	Scale is not level	Level the scale

Fault description	Cause	Remedy
Result is incorrect	Significant variations in ambient temperature	Stabilise the environmental conditions or move the device to another location
	Electromagnetic fields / Electrostatic charge	If possible, switch off the interfering device or change the location of the scale
	Transport safety lock not removed	Remove the transport safety lock ➡ <i>Further information on page 132</i>

18 Product specification

18.1 Technical data

Model-specific data

Model	CDS 4K0.02	CDS 15K0.05	CDS 16K0.1
Item no. / Type	TCDS 4K0.02-A	TCDS 15K0.05-A	TCDS 16K0.1-A
Readout (d)	0,02 g	0,05 g	0,1 g
Weighing capacity (max.)	4 kg	15 kg	16 kg
Reproducibility	0,04 g	0,1 g	0,1 g
Linearity	+/-0,1 g	+/-0,25 g	+/-0,3 g
Stabilization time (typical)	3 s	3 s	3 s
Smallest part weight when piece counting (laboratory conditions)	20 mg	50 mg	100 mg
Smallest part weight when piece counting (normal conditions)	200 mg	500 mg	1 g
Adjustment points	2 / 4 kg	2 / 5 / 10 / 15 kg	2 / 5 / 10 / 15 kg
Recommended adjustment weight (class)	2 kg (F1) / 2 kg (F1)	15 kg (F1)	15 kg (F1)
Time required to reach operating temperature	2 h	2 h	2 h
Measuring system	Strain gauges (DMS)	Strain gauges (DMS)	Strain gauges (DMS)
Device for suspended weighing	No	Hook	Hook
		Optional item (not included in the scope of delivery)	Optional item (not included in the scope of delivery)
Transport safety lock	Yes	No	No

Model	CDS 30K0.1	CDS 30K0.1L	CDS 36K0.2L
Item no. / Type	TCDS 30K0.1-A	TCDS 30K0.1L-B	TCDS 36K0.2L-B

Product specification

Model	CDS 30K0.1	CDS 30K0.1L	CDS 36K0.2L
Readout (d)	0,1 g	0,1 g	0,2 g
Weighing capacity (max.)	30 kg	30 kg	36 kg
Reproducibility	0,2 g	0,2 g	0,2 g
Linearity	+/-0,5 g	+/-0,3 g	+/-0,6 g
Stabilization time (typical)	3 s	3 s	3 s
Smallest part weight when piece counting (laboratory conditions)	100 mg	100 mg	200 mg
Smallest part weight when piece counting (normal conditions)	1 g	1 g	2 g
Adjustment points	10 / 15 / 20 / 30 kg	10 / 20 / 30 kg	12 / 24 / 36 kg
Recommended adjustment weight (class)	30 kg (E2)	30 kg (F1)	36 kg (F1)
Time required to reach operating temperature	2 h	2 h	2 h
Measuring system	Strain gauges (DMS)	Strain gauges (DMS)	Strain gauges (DMS)
Device for suspended weighing	Hook	No	No
	Optional item (not included in the scope of delivery)		
Transport safety lock	No	Yes	Yes

Model	CDS 60K0.2
Item no. / Type	TCDS 60K0.2-B
Readout (d)	0,2 g
Weighing capacity (max.)	60 kg
Reproducibility	0,2 g

Model	CDS 60K0.2
Linearity	+/-0,6 g
Stabilization time (typical)	3 s
Smallest part weight when piece counting (laboratory conditions)	200 mg
Smallest part weight when piece counting (normal conditions)	2 g
Adjustment points	20 / 30 / 50 / 60 kg
Recommended adjustment weight (class)	60 kg (F1)
Time required to reach operating temperature	2 h
Measuring system	Strain gauges (DMS)
Device for suspended weighing	No
Transport safety lock	Yes

Dimensions and weight

Model type (dimensions and weight)	CDS 4K0.02	CDS 15K0.05	CDS 16K0.1
Dimensions of the device (without removable parts)	225 x 115 x 60 mm	225 x 115 x 60 mm	225 x 115 x 60 mm
Platform dimensions	230 x 230 x 103 mm	308 x 318 x 75 mm	308 x 318 x 75 mm
Net weight	8 kg	8 kg	8 kg

Model type (dimensions and weight)	CDS 30K0.1	CDS 30K0.1L	CDS 36K0.2L
Dimensions of the device (without removable parts)	225 x 115 x 60 mm	225 x 115 x 60 mm	225 x 115 x 60 mm
Platform dimensions	308 x 318 x 75 mm	500 x 400 x 127 mm	500 x 400 x 127 mm

Product specification

Model type (dimensions and weight)	CDS 30K0.1	CDS 30K0.1L	CDS 36K0.2L
Net weight	8 kg	12 kg	12 kg

Model type (dimensions and weight)	CDS 60K0.2
Dimensions of the device (without removable parts)	225 x 115 x 60 mm
Platform dimensions	500 x 400 x 127 mm
Net weight	12 kg

General data

Units		kg, g, dwt, tl (TW), tl (HK), ozt, tl (Singap, Malays), ct, mo, lb, oz
Interface		KUP connection (KERN Universal Port) <ul style="list-style-type: none"> ■ RS232 ■ USB-D ■ Ethernet ■ Wi-Fi ■ Bluetooth ■ Analogue
Ambient conditions	Ambient temperature	-10–+40 °C
	Relative humidity (non-condensing)	80 %

Power supply

Scale	Input	Voltage	5,9 V
		Current	1 A
Power supply unit	Input	Voltage	100–240 V AC
		Frequency	50–60 Hz
Rechargeable battery	Charging time		8 h
	Operating temperature range with battery		0–40 °C
	Operating time (without backlight)		48 h
	Operating time (with backlight)		24 h
	Note		Optional item (not included in the scope of delivery)

18.2 Declaration of Conformity

The current EC/EU Declaration of Conformity can be found online at ➔ www.kern-sohn.com/shop/en/downloads.

19 Contact details

KERN Hotlines

Technical questions about our products?

Here you will find quick help: +49 7433 9933- ...

Service hotline

For general technical questions about your KERN product
→ **199**

Laboratory and analytical balances

For all technical questions relating to our high-quality precision balances, analytical balances (especially with force-compensated measuring systems, tuning fork and high application density)
→ **444**

Industrial scales

For all technical questions relating to our basic scales, e.g. pocket balances, school balances, bench scales, platform scales, counting scales, counting systems, floor scales, pallet truck scales, crane scales, veterinary scales
→ **333**

Medical scales

For all technical questions relating to our baby scales, personal floor scales, handrail scales, chair scales, wheelchair scales, hand grip dynamometers
→ **499**

Optical devices

For all technical questions relating to our microscopes, microscope cameras, refractometers
→ **777**

EasyTouch (software)

For all technical questions regarding installation and commissioning of the EasyTouch software
→ **322**

Calibration and verification service

For all questions relating to our calibration and verification services
→ **196**

20 Index

A

Adjusting	141
External adjusting using a specified weight	142
External adjusting using a user-defined weight	143
Application	
Accessing and using the application menu	165
Counting	172
Tolerance weighing	181
Weighing	167
Application menu	
"Counting" overview	172
"Tolerance weighing" overview	181
"Weighing" overview	167
Operating	165
Applications	
Lock applications	147

B

Batteries	
Disposal	192

C

Cleaning	187
Commissioning	138
Contact us	202
Counting	
Check quantity (limit values)	179
Set the reference quantity	174
Simple piece counting	176
Weigh target number of items	177

D

Decommissioning	190
Display	126
Backlighting	146
Disposal	192
Document	113
Other versions	114
Presentation	114

E

Environmental conditions	121
--------------------------	-----

H

Hold function	168
---------------	-----

I

Icons	
Display	126
Install	131
Installation	
Levelling	133
Remove the transport safety lock	132
Safety	131
Select place of use / location	131
Wired interfaces	134
Installation location	121
Interfaces	149
Output settings	151
Settings	150
Setup	149
Wired device connections	134

K

KCP (KERN Communications Protocol)	160
Keys	124
Change key assignments	145

L

Level indicator	133
Levelling	133

M

Mains connection	135
Malfunctions	194

O

Operation	
Multiply the weight value	170
Percentage weighing	169
Simple weighing	167
Switching off	162
Switching on	161
Tare	163
Zero setting	162
Optional software	149
Overview	
"Counting" application menu	172
"Tolerance weighing" application menu	181
"Weighing" application menu	167
Device	123

Display	126	Backlighting	146
Keys	124	Factory settings	148
Output settings	151	Filter	147
Setup menu	140	Interface parameters	150
		Keypad sound	144
		Taring range	146
		Zero Limit	146
		Zero Tracking	147
		Zero-setting range	146
P		Setup menu	
Place of use	121	Overview	140
Power supply unit		Signal words	116
Changing the power plug adapter	134	Storage	190
Connection	135	Symbol	
Pre-tare	164	Warning	117
Presentation conventions	114	T	
Printout		Technical data	197
Output settings	151	Tolerance weighing	
Outputting the totals	186	Check weight (limit values)	184
Set layout	157	Weigh target weight	182
		Totalizing function	186
R		Transport	129
Rechargeable battery	136	Transport safety lock	132
Charge	137	Transport safety lock	132
Insert the rechargeable battery	136	Troubleshooting	194
Safety precautions	118	U	
Restore factory settings	148	Underfloor weighing	167
Return		Units	
Goods return	191	Set available display units	147
Warranty	191	W	
S		Warning signs	117
Safety precautions	116	Warnings	116
Batteries and rechargeable batteries	118		
Improper use	120		
Intended use	120		
IT security	120		
Operators' obligations	120		
Place of use	121		
Power supply units	117		
User qualification	121		
Scope of delivery	129		
Servicing	188		
Maintenance schedule	188		
Settings	139		
Adjusting	141		
Assignment of keys	145		
Automatic shut-off function	144		
Automatic switch-on function	145		
Automatic zero tracking	147		
Available applications	147		
Available display units	147		

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